

CRインレー形成ダイバーセット

医療機器届出番号 09B1X00006002010
一般医療機器 一般の名称：歯科用ダイヤモンドバー
販売名 マニーダイヤモンド









作業手順の短縮

CRインレー窩洞は円形窩洞にする事をお奨めします。

円形窩洞が特長！



サイズ：38(H)×75(W)×60(L)mm
最小包装単位：6本入/セット

咬合面形成用バー					隣接面形成用バー		
	 <p>ISO 196/020 5.0/19.0 CR-11F</p>	 <p>ISO 196/019 5.0/19.0 CR-11EF</p>	 <p>ISO 196/015 4.0/19.0 CR-12F</p>	 <p>ISO 196/014 4.0/19.0 CR-12EF</p>		 <p>ISO 197/012 7.0/21.5 CR-21F</p>	 <p>ISO 289/013 8.0/21.5 CR-22F</p>