

## MICRO FILE Type-K (MICROSCOPES INSTRUMENTS)

---

**MICRO FILE Type K(Taper:2/100)**

**6pcs/case**

**Assortment:15~40**

**Working Length:(L)**

### HOW TO USE

● Fit the file in the ENDO HOLDER (See Below) then use it with filing movement.

## ENDO HOLDER (MICROSCOPES INSTRUMENTS)

---

**Screw-on, Exclusive Holder for MICRO FILES, SPREADERS, DIA FILES.**

### HOW TO USE

● Be fitted with the MICRO FILES, MICRO SPREADERS, DIA