

ULTRASONIC ENDO FILE (MICROSCOPES INSTRUMENTS)

MICRO ENDO FILE(Taper:5/100)

6pcs/case

Working Length:(L)

HOW TO USE

- **Fit the file in the ENDO HOLDER (See Below) then use it with an ultrasonic scaler.**
- **Scaler power should be less than the rate "3" for your**

ENDO HOLDER (MICROSCOPES INSTRUMENTS)

Screw-on, Exclusive Holder for MICRO FILES, SPREADERS, DIA FILES.

HOW TO USE

- **Be fitted with the MICRO FILES, MICRO SPREADERS, DIA**