

カーバイドバー

医療機器届出番号 09B1X00006003010
 一般医療機器 一般的名称：歯科用カーバイドバー
 販売名 マニーカーバイドバー

特徴


●刃部とシャンク部の一体形成

マニーカーバイドバーは、刃部とシャンク部の一体形成により、切削中の破折を最小限に抑えました。（＃6・7・8・703・7408は、接合タイプ）

●主な使用


■ROUND

齶窩の開拡・形成に有効



■INVERTED CONE


アンダーカット形成に適しています




































※クロスカットが多いほど
 切削効率が高まります。










●材質 超硬、ステンレススチール










■形態表の見方
















ROUND (PLAIN)							
 0.50 19.30	 0.60 19.30	 0.80 19.30	 1.00 19.30	 1.20 19.30	 1.40 19.30	 1.60 19.30	 1.80 19.30
1/4	1/2	1	2	3	4	5	6
001/005	001/006	001/008	001/010	001/012	001/014	001/016	001/018
INVERTED CONE (PLAIN)		INVERTED CONE (PLAIN)				CYLINDER FLAT END (PLAIN)	
 2.10 19.30	 2.30 19.30	 0.70 19.30	 0.80 19.30	 1.10 19.30	 1.20 19.30	 1.40 19.30	 4.00 19.30
7	8	33 1/2	34	35	36	37	56
001/020	001/022	010/006	010/008	010/010	010/012	010/014	107/009
TAPER FLAT END (PLAIN)				PEAR (PLAIN)			
 4.00 19.30	 4.00 19.30	 5.00 21.00	 4.00 19.30	 4.00 19.30	 4.00 19.30	 3.00 19.30	 4.00 19.30
57	58	169L	170	171	172	245	246
107/010	107/012	171/009	168/010	168/012	168/016	237/009	234/011

TAPER FLAT END(CROSSCUT)								
								
2.00 19.30	2.00 19.30	2.00 19.30	4.00 19.30	4.00 19.30	4.00 19.30	4.00 19.30	4.00 19.30	4.00 19.30
330	331	332	556	557	558	559	560	560
232/008	232/010	232/012	107/009	107/010	107/012	107/014	107/016	107/016

TAPER FLAT END(CROSSCUT)					CYLINDER ROUND END(PLAIN)			
								
4.00 19.30	4.00 19.30	4.00 19.30	7.00 21.00	4.00 19.30	4.70 19.30	4.00 19.30	4.00 19.30	4.00 19.30
699	700	701	701L	702	703	1156	1157	1157
168/009	168/010	168/012	171/012	168/016	168/021	137/009	137/010	137/010

CYLINDER ROUND END(CROSSCUT)					SPECIAL SHAPE			
								
4.00 19.30	4.00 19.30	4.00 19.30	4.00 19.30	4.00 19.30	5.50 21.00	6.00 21.00	8.9 23.4	8.9 23.4
1158	1159	1557	1558	1559	M16DB	M14CF	CE4CF	CE4CF
137/012	137/014	137/010	137/012	137/014	196/016	164/014	166/014	166/014

FINISHING(PLAIN)							
							
4.00 19.30	4.00 19.30	4.00 19.30	9.00 21.00	11.00 23.00	4.00 20.00	4.00 19.30	4.00 19.30
701X	M21DC	M23DC	M15EZ	M15ZB	7408	7901	7902
168/012	137/010	168/012	198/015	199/015	277/023	496/009	496/010

REMOVAL(CROSSCUT)				
				
4.00 19.30	2.00 19.00	4.00 19.00	4.00 19.00	4.00 19.00
7903	1931	1957	1958	1970
496/012	237/010	137/010	137/012	194/010

● 最小包装単位 ○Ⅰクラス 10本入/シート
 △Ⅱクラス 10本入/シート
 ◎Ⅲクラス 5本入/シート

■ 標準価格 (1シート)

○ Ⅰクラス △ Ⅱクラス
¥2,200 **¥2,800**

◎ Ⅲクラス
¥4,600