



Healthy teeth, happy smiles,  
to the world.

**MANI**



---

**D I A - B U R S**  
**Endo Access Burs**  
**S E L E C T I O N**

---

Basic and essential burs for root canal access



	Pulp chamber and Chamber wall		Internal axial wall	Canal openings
Item code	EA-10	EA-L10	EA-13	EA-40
Grid size	Standard		Standard	Standard
Total Length (mm)	22.0	25.0	28.0	25.0
Working Length (mm)	11.0	11.0	10.2	1.4
Max. Diameter (mm)	1.9	1.9	1.8	1.6
Max. Allowable Speed (K)	300	300	90	160
Packaging	Strip of 5		Strip of 5	Strip of 5