







Symbols Glossary

Symbol	Symbol Title	Description of Symbol	Standard Reference
	Manufacturer	Indicates the medical device manufacturer	ISO 15223-1:2016 5.1.1 ISO 15223-1:2021 5.1.1
	Authorized representative in the European Community/ European Union	Indicates the authorized representative in the European Community/ European Union	ISO 15223-1:2016 5.1.2 ISO 15223-1:2021 5.1.2
	Date of manufacture	Indicates the date when the medical device was manufactured	ISO 15223-1:2016 5.1.3 ISO 15223-1:2021 5.1.3
	Use-by date	Indicates the date after which the medical device is not to be used	ISO 15223-1:2016 5.1.4 ISO 15223-1:2021 5.1.4
	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified	ISO 15223-1:2016 5.1.5 ISO 15223-1:2021 5.1.5
	Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified	ISO 15223-1:2016 5.1.6 ISO 15223-1:2021 5.1.6




Symbols Glossary

Symbol	Symbol Title	Description of Symbol	Standard Reference
	Importer	Indicates the entity importing the medical device into the locale	ISO 15223-1:2021 5.1.8
	Sterilized using ethylene oxide	Indicates a medical device that has been sterilized using ethylene oxide	ISO 15223-1:2016 5.2.3 ISO 15223-1:2021 5.2.3
	Do not resterilize	Indicates a medical device that is not to be resterilized	ISO 15223-1:2016 5.2.6 ISO 15223-1:2021 5.2.6
	Non-sterile	Indicates a medical device that has not been subjected to a sterilization process	ISO 15223-1:2016 5.2.7 ISO 15223-1:2021 5.2.7
	Do not use if package is damaged and consult instructions for use	Indicates that a medical device that should not be used if the package has been damaged or opened and that the user should consult the instructions for use for additional information	ISO 15223-1:2016 5.2.8 ISO 15223-1:2021 5.2.8
	Single sterile barrier system	Indicates a single sterile barrier system	ISO 15223-1:2021 5.2.11

Symbols Glossary

Symbol	Symbol Title	Description of Symbol	Standard Reference
	Do not re-use	Indicates a medical device that is intended for one single use only	ISO 15223-1:2016 5.4.2 ISO 15223-1:2021 5.4.2
	Consult instructions for use or consult electronic instructions for use	Indicates the need for the user to consult the instructions for use	ISO 15223-1:2016 5.4.3 ISO 15223-1:2021 5.4.3
	Contains hazardous substances	Indicates a medical device that contains substances that can be carcinogenic, mutagenic, reprotoxic (CMR), or substances with endocrine disrupting properties	ISO 15223-1:2016 5.4.10 ISO 15223-1:2021 5.4.10
	Medical device	Indicates the item is a medical device	ISO 15223-1:2021 5.7.7
	Unique device identifier	Indicates a carrier that contains unique device identifier information	ISO 15223-1:2021 5.7.10
Rx Only Or 	Prescription use only	Caution: Federal law restricts this device to sale by or on the order of a physician	21 CFR 801.109

Symbols Glossary

Symbol	Symbol Title	Description of Symbol	Standard Reference
	CE mark	European conformity	The requirements for accreditation and market surveillance relating to the marketing of products; Medical Device Directive & Medical Device Regulation.
 Or 	CE mark with Notified Body Number	European conformity	The requirements for accreditation and market surveillance relating to the marketing of products; Medical Device Directive & Medical Device Regulation.